

# Application Report 11

## Industrial Solvents on SUPELCOWAX 10

In this application we analyzed a mix of 58 commonly used industrial solvents on the SUPELCOWAX 10. This polar column offered excellent peak shape, response and resolution and would be beneficial to use in combination with a nonpolar column such as the Equity-1 for this application.


### Key Words

solvents, SUPELCOWAX, 24211, polar

Author: K. Stenerson

Raw Data File Name:  
GS100APPSWAX0213.D

Acquisition System: GC 6248



### Conditions

Column: SUPELCOWAX 10, 30m x 0.32mm ID, 1.0µm  
Cat. No.: 24211  
Oven: 35°C (8 min) to 130°C @ 4°C/min (2 min)  
Inj.: 250°C  
Det.: FID, 250°C  
Flow: Helium, 25cm/sec constant @ 35°C  
Injection: 0.5µL, split 200:1  
Liner: Split cup design  
Sample: 0.5µL of a 59 component neat solvent mixture

### Peak IDs

- |                                  |                                  |
|----------------------------------|----------------------------------|
| 1. Methanol                      | 31. n-Propyl acetate             |
| 2. Methyl formate                | 32. 4-Methyl-2-pentanone         |
| 3. Ethanol                       | 33. Isoamyl alcohol              |
| 4. Acetone                       | 34. Dimethylformamide            |
| 5. 2-Propanol                    | 35. Toluene                      |
| 6. Ethyl formate                 | 36. Isobutyl acetate             |
| 7. 1,1-Dichloroethylene          | 37. 2-Hexanone                   |
| 8. Methylene chloride            | 38. Mesityl oxide                |
| 9. Methyl acetate                | 39. Tetrachloroethene            |
| 10. 1-Propanol                   | 40. n-Butyl acetate              |
| 11. trans-1,2-Dichloroethylene   | 41. Diacetone alcohol            |
| 12. 1,1-Dichloroethane           | 42. Chlorobenzene                |
| 13. 2-Butanone                   | 43. 5-Methyl-2-hexanone          |
| 14. sec-Butanol                  | 44. Ethyl benzene                |
| 15. Hexane                       | 45. m-Xylene                     |
| 16. Ethyl acetate                | 46. p-Xylene                     |
| 17. Chloroform                   | 47. Isoamyl acetate              |
| 18. Tetrahydrofuran              | 48. Cyclohexanol                 |
| 19. Isobutanol                   | 49. Styrene                      |
| 20. 2-Methoxyethanol             | 50. o-Xylene                     |
| 21. 1,2-Dichloroethane           | 51. 1,1,2,2-Tetrachloroethane    |
| 22. 1,1,1-Trichloroethane        | 52. 2-Ethoxyethyl acetate        |
| 23. Isopropyl acetate            | 53. Butyl cellosolve             |
| 24. n-Butanol                    | 54. n-Amyl acetate isomers       |
| 25. Benzene                      | 55. 2-Methylcyclohexanol isomers |
| 26. Carbon tetrachloride         | 56. 1,2-Dichlorobenzene          |
| 27. 2-Nitropropane               | 57. 2-Methylphenol               |
| 28. Trichloroethylene            | 58. 3-Methylphenol               |
| 29. 1,4-Dioxane                  | 59. 4-Methylphenol               |
| 30. 2-Ethoxyethanol (cellosolve) |                                  |