

Technical Bulletin

Procainamide Labeled Dextran Ladder

Catalogue number SMB01378

Product Description

Glycans are associated with numerous biological processes due to their structural and binding properties. New functional roles of glycans are constantly being discovered. However, many glycans are not assigned with any functional aspect. To better understand their roles, glycan profiling is essential.

Procainamide labelling of glycans is one method of profiling glycans by coupling with hydrophilic interaction liquid chromatography (HILIC) and fluorescence detection (FD) or mass spectrometry (MS).

The procainamide labeled dextran provides for assignment of GU retention values of unknown glycans to aid in identification and standardization over time.

Key benefits are listed below:

1. Each labeled glycan - Glucose Unit (GU) value increases in an incremental value of its GU thus helping an easy assignment of unknown glycans.
2. Assignment of sample glycans can be done against dextran ladder. GU retention values allows different instruments to be used, and sample comparisons can be made over time for quality control.
3. The observed GU values are between 2 to 30, thus enabling larger glycans retention profiling.
4. The ladder can also serve as a system suitability standard for each batch of samples.
5. Elution of glycans in ladder happens with increasing units, thus enabling the use on day-to-day analysis in LC.
6. A polynomial fit can then be used to assign RT values of sample glycans for identification or confirmation in a QC setting.

Components of procainamide labeled dextran ladder

*Procainamide labeled dextran ladder containing increasing GU units. Each vial contains 200 µg of lyophilized powder.

Precautions and Disclaimer

This product is for R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Store the sealed vials at 4 °C.

Preparation Instructions

The lyophilized material consists of a broad range of both straight chains, α -1,6 glucose linkages, and larger branched chain oligosaccharides with α -1,3 linkages. The use of 75% acetonitrile (water) solubilizes the shorter chain GU components while use of 50% acetonitrile in water will solubilize many of the longer, branched chain components as well. Depending on the GU range of interest one may choose more, or less, acetonitrile for the reconstitution. See Figure 1 below.

One vial of procainamide labeled dextran ladder contains 200 µg of lyophilized powder. For reconstitution, dissolve the sample in 200 µL of 75% acetonitrile. Mix well and centrifuge to remove any undissolved material.

Note: Dilution can be adjusted depending on the needs of experiments.

Procedure

LC Conditions:

Column: BIOshell™ Glycan HPLC Column
(150 mm x 2.1 mm x 2.7 µm) 90 Å.

Column oven temperature: 40 °C

Sample temperature: 5 °C

Flow rate: 0.3 mL/min

Eluent A: 75 mM Ammonium formate in water

Eluent B: ACN

Injection volume: 10 µL

Fluorescence detector: Ex: 308 nm Em: 359 nm

Gradient:

Time (min)	A [%]	B [%]
0	20	75
75	41	59


Figure 1: Chromatogram showing procainamide labeled dextran ladder, as measured by fluorescence. The solvent used to reconstitute the sample has an influence on the series of dextran components seen.

For details on labelling of released glycans using procainamide in the laboratory, please see the application note:

[Step-by-Step Protocol for Procainamide Labelled Glycan Profiling of a Monoclonal Antibody](#)

References

1. Klapoetke, S., *et al.*, The evaluation of a novel approach for the profiling and identification of N-linked glycan with a procainamide tag by HPLC with fluorescent and mass spectrometric detection. *J Pharm Biomed Anal*, **53(3)**, 315-24 (2010).
2. HPLC Analysis of a Procainamide-Labeled Dextran Ladder on BIOshell™ Glycan using HILIC-FLR.
(<https://www.sigmaaldrich.com/IN/en/technical-documents/chromatograms/hplc/hplc-analysis-of-a-procainamide-labeled-dextran-ladder-on-bioshell-glycan-using-hilic-flr/supelco/q006408>)



Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

smb01378 Rev 09/24

**MILLIPORE
SIGMA**