



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

4-Methylumbelliferyl p-trimethylammonio cinnamate chloride

Product Number **M 4507**
Storage Temperature -0 °C

Product Description

Molecular Formula: $C_{22}H_{22}ClNO_4$
Molecular Weight: 399.9
CAS Number: 34197-47-2
Melting Point: 204-206 °C¹
Synonym: MUTMAC

4-Methylumbelliferyl p-trimethylammoniocinnamate chloride (MUTMAC) is a good titrant for α -chymotrypsin. When α -chymotrypsin reacted with MUTMAC, acylation was complete in less than 2 min at pH values between 7.4 and 8.3, and the rate of deacylation was negligible. A solution of α -chymotrypsin was assayed spectrofluorimetrically with MUTMAC. The released 4-methylumbelliferone was measured using excitation and emission wavelengths of 365 nm and 410 nm, respectively.¹ The amount of 4-methylumbelliferone liberated was proportional to the enzyme concentration (12 μ M) and independent of both the concentration of MUTMAC (0.16 to 1.600 μ M) and of pH in the range 6.8 to 8.3. The MUTMAC reacts with α -chymotrypsin in less than 2 minutes, has no turnover and is sensitive to 10^{-11} mole of enzyme.²

Solutions of α - and β -trypsin, thrombin, Factor X_a , α -chymotrypsin can be titrated spectrofluorimetrically with MUTMAC with accuracy and precision comparable to spectrophotometric methods. The spectrofluorimetric method can be used with as little as 0.02 nmol of enzyme; the spectrophotometric method requires about 1 nmol of enzyme.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is soluble in pyridine:water (1:1, 50 mg/ml) with heat, yielding a clear, light yellow/green solution. MUTMAC has also been dissolved in water (40 mg/ml).¹

Storage/Stability

Aqueous solutions (0.05-0.2 mM, 20-80 μ g/ml) of MUTMAC can be stored frozen for several weeks.¹

References

1. Jameson, G. W., et al., Determination of the operational molarity of solutions of bovine α -chymotrypsin, trypsin, thrombin and factor X_a by spectrofluorimetric titration. *Biochem. J.*, **131**(1), 107-117 (1973).
2. Coleman, P. L., et al., Some sensitive methods for the assay of trypsinlike enzymes. *Meth. Enzymol.*, **45**, 12-26 (1976).

MES/CRF 1/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.