

Technical Bulletin

Lectin from *Phaseolus vulgaris* (red kidney bean)

Phytohemagglutinin PHA-P, lyophilized powder

L8754

Product Description

Lectins are proteins or glycoproteins of non-immune origin that agglutinate cells and/or precipitate complex carbohydrates. Lectins can bind glycoproteins even in the presence of various detergents.¹ The agglutination activity of these highly specific carbohydrate-binding molecules is usually inhibited by a simple monosaccharide. However, for some lectins, di-, tri-, and even polysaccharides are required.

Lectins are isolated from a wide variety of natural sources, including seeds, plant roots and bark, fungi, bacteria, seaweed and sponges, mollusks, fish eggs, body fluids of invertebrates and lower vertebrates, and from mammalian cell membranes.² Lectins have proven very valuable in a wide variety of applications *in vitro*, including:

- blood grouping and erythrocyte polyagglutination studies
- mitogenic stimulation of lymphocytes
- lymphocyte subpopulation studies
- fractionation of cells and other particles
- histochemical studies³

A range of lectins suitable for the above applications is available. Most Sigma-Aldrich® lectins are highly purified by affinity chromatography. However, some are offered as purified or partially purified lectins, suitable for specific applications. Many of the lectins are available conjugated to the following moieties, where conjugation does not alter the specificity of the lectin:

- fluorochromes (for detection by fluorimetry)
- enzymes (for enzyme-linked assays)
- insoluble matrices (for use as affinity media)

Please refer to Table 1 for general information on the most common lectins.

Lectin PHA-P is not inhibited easily by monosaccharides, but may be inhibited by oligosaccharides. PHA-P is a mixture of PHA-E (MW = 128 kDa) and PHA-L (MW = 126 kDa).

Several publications,^{4,5} theses^{6,7} and dissertations⁸⁻¹¹ cite use of product L8754 in their research protocols.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Aggregation is thought to occur in the presence of high concentrations of 2-mercaptoethanol.

Preparation Instructions

This lectin is soluble in phosphate buffered saline, pH 6.8, at 1 mg/mL.

Procedure

A general agglutination procedure using this lectin with 96-well plates is as follows:

1. Prepare a lectin solution of 1 mg/mL in PBS buffer, pH 6.8.
2. Pipette 50 µL of fresh PBS into each well and add 50 µL of the lectin solution into the first well.
3. Serial dilutions are made by pipetting 50 µL from each successive well into the next well.
4. Blood type A with a 2% hematocrit is used as the substrate.
5. Pipette 50 µL of blood into each well.
6. Visually determine agglutination.

References

1. Lotan, R. *et al.*, *Biochemistry*, **16(9)**, 1787-1794 (1977).
2. Chrispeels, M.J., and Raikhel, N.V., *Plant Cell*, **3(1)**, 1-9 (1991).
3. Schumacher, U. *et al.*, "Lectins as tools in histochemical techniques: a review of methodological aspects", in *Lectin Reviews* (New York), Vol. 1, pp. 195-201 (1991).
4. Plasman, M. *et al.*, *Biol. J. Linnean Soc.*, **128(4)**, 936-943 (2019).
5. Rekdal, S.L. *et al.*, *J. Avian Biol.*, **52(5)**, doi.org/10.1111/jav.02734 (2021).
6. Dorton, Kristy L., "Effects of copper source and concentration on copper status and immune function in growing and finishing steers". Colorado State University, M.S. thesis, p. 38 (2002).
7. Lazzaroni, Sharon Marie, "Immune responses in *Burkholderia pseudomallei* infection, relapse, reactivation and protection". James Cook University, M.Sc. thesis, p. 131 (2009).
8. Yang, David Shih-Fan, "Immunomodulatory effects of interferon-alpha on T cell subsets in Behcet's disease *ex vivo*, and the *in vitro* effects of treatment on healthy donor cells". University College London, Ph.D. dissertation, p. 113 (2011).
9. Rive, Craig Michael, "The Immunopathogenesis of drug hypersensitivity". Murdoch University, Ph.D. dissertation, pp. 37, 47 (2016).
10. Zeytin, Sinem, "Towards optimisation of fish larval nutrition: the effects of new components in microdiet composition and feeding practices on larval digestive physiology". Christian-Albrechts-Universität Kiel, Dr. rer. nat. dissertation, p. 34 (2018).
11. Rieger, Jan Hendrik, "Targets of tumor-infiltrating lymphocyte reactivity in pancreatic cancer". Ruperto Carola University Heidelberg, Dr. rer. nat. dissertation, p. 28 (2019).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

Table 1. General lectins guide

Lectin	MW (kDa)	Subunits	Specificity		Mitogenic Activity
			Blood Group	Sugar	
<i>Abrus precatorius</i>					+
• Agglutinin	134	4		Gal	
• Abrin A (toxin)	60	2		Gal	
• Abrin B (toxin)	63.8	2($\alpha\beta$)		Gal	
<i>Agarius bisporus</i>	58.5	-	-	β -gal(1 \rightarrow 3)galNAc	+
<i>Anguilla anguilla</i>	40	2	H	α -L-Fuc	
<i>Arachis hypogaea</i>	120	4	T	β -gal(1 \rightarrow 3)galNAc	
<i>Artocarpus integrifolia</i>	42	4	T	α -gal \rightarrow OMe	+
<i>Bandeiraea simplicifolia</i>					
• BS-I	114	4	A, B	α -gal, α -galNAc	
• BS-I-A ₄	114	4	A	α -galNAc	
• BS-I-B ₄	114	4	B	α -gal	
• BS-II	113	4	acq, B, Tk, T	glcNAc	
<i>Bauhinia purpurea</i>	195	4		β -Gal(1 \rightarrow 4)GalNAc	+
<i>Caragana arborescens</i>	60; 120 ^a	2/4	-	GalNAc	
<i>Cicer arietinum</i>	44	2	-	Fetuin	
<i>Codium fragile</i>	60	4	-	GalNAc	
Concanavalin A	102	4	-	α -Man, α -Glc	+
Succinyl-Concanavalin A	51	2	-	α -Man, α -Glc	+ ^b
<i>Cytisus scoparius</i>	-	-	-	GalNAc, Gal	
<i>Datura stramonium</i>	86	2($\alpha\beta$)	-	(GlcNAc) ₂	
<i>Dolichos biflorus</i>	140	4	A ₁	α -GalNAc	
<i>Erythrina corallodendron</i>	60	2	-	β -Gal(1 \rightarrow 4)GlcNAc	+
<i>Erythrina cristagalli</i>	56.8	2($\alpha\beta$)	-	β -Gal(1 \rightarrow 4)GlcNAc	
<i>Euonymus europaeus</i>	166	4($\alpha\beta$)	B,H	α -Gal(1 \rightarrow 3)Gal	+
<i>Galanthus nivalis</i>	52	4	(h)	Non-reduced α -Man	
<i>Glycine max</i>	110	4	-	GalNAc	+ ^c
<i>Helix aspersa</i>	79	-	A	GalNAc	
<i>Helix pomatia</i>	79	6	A	GalNAc	
<i>Lathyrus odoratus</i>	40-43	4($\alpha\beta$)	-	α -Man	+
<i>Lens culinaris</i>	49	2	-	α -Man	+
<i>Limulus polyphemus</i>	400	18	-	NeuNAc	
• Bacterial agglutinin	-	-	-	GalNAc, GlcNAc	
<i>Lycopersicon esculentum</i>	71	-	-	(GlcNAc) ₃	
<i>Maackia amurensis</i>	130	2($\alpha\beta$)	O	Sialic acid	+
<i>Maclura pomifera</i>	40-43	2($\alpha\beta$)		α -Gal, α -GalNAc	
<i>Momordica charantia</i>	115-129	4($\alpha\beta$)	-	Gal, GalNAc	
<i>Naja mocambique mocambique</i>	-	-	-	-	
<i>Naja naja kaouthia</i>	-	-	-	-	

<i>Narcissus pseudonarcissus</i>	26	2	(h)	α -D-Man	
<i>Persea americana</i>	-	-	-	-	
<i>Phaseolus coccineus</i>	112	4	-	-	
<i>Phaseolus limensis</i>	247 (II)	8	A	GalNAc	+
	124 (III)	4			
<i>Phaseolus vulgaris</i>					
PHA-E	128	4	-	Oligosaccharide	+
PHA-L	128	4	-	Oligosaccharide	+
PHA-P					
PHA-M					
<i>Phytolacca americana</i>	32	-	-	(GlcNAc) ₃	+
<i>Pisum sativum</i>	49	4($\alpha\beta$)	-	α -Man	+
<i>Pseudomonas aeruginosa PA-I</i>	13 – 13.7	-	-	Gal	+ ^c
<i>Psophocarpus tetragonolobus</i>	35	1	-	GalNAc, Gal	
<i>Ptilota plumosa</i>	65; 170	-	B	α -Gal	
<i>Ricinus communis</i>					
• Toxin, RCA60	60	2	-	GalNAc, β -Gal	
• Toxin, RCA120	120	4	-	β -Gal	
<i>Sambucus nigra</i>	140	4($\alpha\beta$)	-	α NeuNAC(2→6)Gal, GalNAc	+ ^c
<i>Solanum tuberosum</i>	50, 100 ^a	1,2	-	(GlcNAc) ₃	
<i>Sophora japonica</i>	133	4	A,B	β -GalNAc	
<i>Tetragonolobus purpureas</i>	120 [A]	4	H	α -L-Fuc	
	58 [BA]	2	H	α -L-Fuc	
	117 [C]	4	H	α -L-Fuc	
<i>Triticum vulgare</i>	36	2	-	(GlcNAc) ₂ , NeuNAC	+
<i>Ulex europaeus</i>					
UEA I	68	-	H	α -L-Fuc	
UEA II	68	-		(GlcNAc) ₂	
<i>Vicia faba</i>	50	4($\alpha\beta$)	-	Man,Glc	+
<i>Vicia sativa</i>	40	4($\alpha\beta$)	-	Glc,Man	+
<i>Vicia villosa</i>	139	4	A ₁ +T _n	GalNAc	
• A4	134	4	A ₁	GalNAc	
• B4	143	4	T _n	GalNAc	
<i>Vigna radiata</i>	160	4	-	α -Gal	
<i>Viscum album</i>	115	4($\alpha\beta$)	-	β -Gal	
<i>Wisteria floribunda</i>	68	2	-	GalNAc	

^a Concentration-dependent molecular weight

^c Mitogenic for neuraminidase-treated lymphocyte

^b Non-agglutinating and mitogenic

The life science business of Merck operates
as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates.
All other trademarks are the property of their respective owners. Detailed information on
trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

L8754pis Rev 05/22 IRB,MWM,JRC,NSB,SAG,MAP,GCY

