

## Product Information

# Streptavidin-Agarose from *Streptomyces avidinii*

Buffered aqueous suspension

**S1638**

## Product Description

Streptavidin is a 66 kDa homotetrameric protein, isolated from *Streptomyces avidinii*, which, like avidin, has a high affinity for biotin ( $K_a \sim 10^{15} \text{ M}^{-1}$ ).<sup>1-4</sup> Streptavidin is slightly anionic (pI  $\sim$  5-6) and non-glycosylated. These properties contribute to its relatively low non-specific binding compared to egg white avidin (a glycoprotein with pI  $\sim$  10.5).<sup>1-3,5</sup> Streptavidin is also more resistant than avidin to dissociation into subunits by guanidinium chloride.<sup>6</sup>

Streptavidin-agarose can be used to immobilize or isolate various biotinylated macromolecules and complexes (proteins,<sup>7-9</sup> antibodies,<sup>3,9-12</sup> lectins,<sup>9,13</sup> nucleic acids,<sup>9,14-16</sup> receptors and ligands<sup>9,14</sup>). The inherent high affinity streptavidin-biotin interaction requires harsh conditions to release biotinylated macromolecules. This feature makes streptavidin-agarose useful in a variety of affinity purification applications.<sup>9,14</sup>

This streptavidin-agarose conjugate is prepared by attaching streptavidin to a CNBr-activated 4% agarose through an amino group, to give a 7-atom spacer. Several publications<sup>17-29</sup> and dissertations<sup>30-39</sup> have cited use of S1638 in their research protocols.

## Product Profile

This product is provided as a suspension in 0.01 M sodium phosphate, pH 7.2, containing 0.15 M NaCl and 0.02% sodium azide.

## Storage/Stability

For long-term storage, store the product at 2-8 °C. **DO NOT FREEZE.** Freezing will damage the bead structure.

## Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

## Preparation Instructions

The agarose beads must be washed thoroughly with 5-10 column volumes of buffer, to remove the sodium azide in the storage buffer.

## Sample Procedures

**Note:** These procedures are provided as general guidelines. The end user should empirically determine appropriate and optimal conditions for the specific application.

### Procedure for Column Purification of Biotinylated Proteins

1. Pour the streptavidin-agarose into an appropriate column. Wash with 5-10 column volumes of PBS.
2. Apply the sample that contains the biotinylated protein.
3. Wash with PBS until the absorbance at 280 nm is minimal.
4. Elute biotinylated proteins either by using 6 M guanidine HCl (pH 1.5-2), or by boiling in 2% SDS with 0.4 M urea.<sup>9</sup>
5. Immediately dialyze or desalt eluted samples, if desired.

### Procedure for Immunoaffinity Column Purification of a Protein<sup>3</sup>

1. Pour the streptavidin-agarose into an appropriate column. Wash with 5-10 column volumes of PBS.
2. Apply the biotinylated antibody.
3. Wash the column with PBS until the absorbance at 280 nm is  $<0.01$ - $0.02$ .
4. Apply the sample (antigen) to the column.
5. Wash with PBS until the absorbance at 280 nm is minimal.
6. Elute the sample (antigen) with 0.1 M acetic acid or 0.1 M glycine HCl (pH 2.5), or other elution buffer to dissociate the antibody-antigen interaction.<sup>3,9,16,40</sup> Immediately neutralize the eluted samples.

**Notes:**Removal of biotinylated substances from streptavidin-agarose:

- To avoid the harsh and denaturing conditions required to dissociate the biotin-streptavidin interaction, an alternative option is the use of a monomeric avidin-agarose, which requires milder dissociation conditions, such as the use of free biotin or low pH.<sup>41</sup>
- Another option is to iminobiotinylate the protein of interest. Iminobiotin has a pH-dependent interaction with streptavidin which allows tight binding at pH > 9.5 and dissociation at pH 4.<sup>42,43</sup>
- For biotinylated nucleic acid probes, DTT or 2-mercaptoethanol can reportedly be used to disrupt the streptavidin-biotin interaction.<sup>15</sup>

Binding inhibitors:

Non-dialyzed nonfat dried milk, which may contain biotin<sup>44</sup> and various sugars such as mannose,<sup>45</sup> is reported to interfere with biotin binding to streptavidin and avidin.

Reduction of non-specific binding:

If non-specific binding is a problem, consider the following options:

- Use a different buffer system, or increase the concentration of NaCl to 0.5 M.<sup>46</sup>
- Increase detergent concentration or use a different detergent.
- Use a cleared cell lysate, by centrifuging the cell lysate to remove unwanted components.<sup>16</sup>
- Wash the resin to remove non-specifically bound proteins, such as with a high pH buffer wash (borate buffer, pH 8-9)<sup>46</sup> followed with a low pH buffer wash (acetate buffer, pH 4).
- Consider using a product with a different charge, such as avidin-agarose or monomeric avidin-agarose.

**References**

- Chaiet, I., and Wolf, F.J., *Arch. Biochem. Biophys.*, **106**, 1-5 (1964).
- Bayer, E.A. et al., *Methods Enzymol.*, **184**, 80-89 (1990).
- Bayer, E.A., and Wilchek, M., *Methods in Molecular Biology • 10: Immunochemical Protocols* (M.M. Manson, ed.). Humana Press (Totowa, NJ), Chapters 13-15, pp. 137-162 (1992).
- Argaraña, C.E. et al., *Nucleic Acids Res.*, **14(4)**, 1871-1882 (1986).
- Green, N.M., *Methods Enzymol.*, **184**, 51-67 (1990).
- Green, N.M., *Adv. Protein Chem.*, **29**, 85 (1975).
- Swack, J.A. et al., *Anal. Biochem.*, **87(1)**, 114-126 (1978).
- von Boxberg, Y. et al., *Eur. J. Biochem.*, **190(2)**, 249-256 (1990).
- Wilchek, M. and Bayer, E.A., *Methods Enzymol.*, **184**, 14-45 (1990).
- Updyke, T.V., and Nicolson, G.L., *Methods Enzymol.*, **121**, 717-725 (1986).
- Updyke, T.V. and Nicolson, G.L., *J. Immunol. Methods*, **73(1)**, 83-95 (1984).
- Gretch, D.R. and Stinski, M.F., *Anal. Biochem.*, **163(1)**, 270-277 (1987).
- Buckie, J.W., and Cook, G.M., *Anal. Biochem.*, **156(2)**, 463-472 (1986).
- Diamandis, E.P., and Christopoulos, T.K., *Clin. Chem.*, **37(5)**, 625-636 (1991).
- Jenne, A., and Famulok, M., *BioTechniques*, **26(2)**, 249-254 (1999).
- Ausubel, F.M. et al. (eds.), *Short Protocols in Molecular Biology* (4th edition). John Wiley & Sons, Inc., pp. 10-89, 12-28, 17-5, (1999).
- Lau, P.N.I., and Cheung, P., *Nucleic Acids Res.*, **41(3)**, e49 (2013).
- Liu, C. et al., *Am. J. Physiol. Gastrointest. Liver Physiol.*, **307(7)**, G749-G759 (2014).
- Tu, K. et al., *Hepatology*, **61(1)**, 361-374 (2015).
- Perrody, E. et al., *eLife*, **5**, e19083 (2016).
- Wang, J. et al., *Autophagy*, **12(10)**, 1931-1944 (2016).
- Wallot-Hieke, N. et al., *Autophagy*, **14(5)**, 743-763 (2018).
- Wang, J. et al., *Proc. Nat. Acad. Sci. USA*, **115(20)**, E4604-E4612 (2018).
- Fletcher-Jones, A. et al., *eLife*, **8**, e44252 (2019).
- Tong, Z. et al., *J. Biol. Chem.*, **294(51)**, 19667-19682 (2019).
- Guo, Q.Q. et al., *Front. Mol. Biosci.*, **7**, 609828 (2021).
- Wang, X. et al., *Cell Metab.*, **33(2)**, 350-366.e7 (2021).
- Shmakova, A. et al., *Cells*, **10(6)**, 1425 (2021).
- Khudaberdiev, S. et al., *RNA Biol.*, **18(9)**, 1252-1264 (2021).
- Schriemer, David Christopher, "Development and Application of Laser-based Ionization Techniques in Time-of-flight Mass Spectrometry". University of Alberta, Ph.D. dissertation, p. 160 (1997).

31. Chrostek-Grashoff, Anna Magdalena, "The role of Rac1 in the epidermis and in the hair follicle". Ludwigs-Maximilians-Universität München, Dr. rer. nat. dissertation, p. 76 (2007).
32. Boladeras, Natàlia Grañé, "Novel approaches to understanding hENT2 and hENT2-related proteins: From novel nuclear variants to global networks". Universitat de Barcelona, Ph.D. dissertation, p. 307 (2012).
33. Lau, Nga Ieng (Priscilla), "Histone crosstalks involving H3 phosphorylation and their role in transcriptional regulation". University of Toronto, Ph.D. dissertation, p. 100 (2013).
34. Demagny, Hadrien, "Convergence des Voies de Signalisation Wnt, FGF et TGF-beta au Niveau des Facteurs de Transcription Smad1 et Smad4" ("Convergence of the Wnt, FGF and TGF-beta Signaling Pathways at the Level of the Smad1 and Smad4 Transcription Factors"). Université Pierre et Marie Curie, Paris VI, Ph.D. dissertation, p. 112 (2014).
35. Ng, Fu Liang, "The Impact of the Blood Pressure-Associated Genetic Locus at SLC4A7 on Gene Expression and Intracellular pH Regulation". Queen Mary, University of London, Ph.D. dissertation, p. 124 (September 2016).
36. Brauns-Schubert, Priska, "Investigation of the Regulation of the Histone Acetyltransferase Tip60 by Phosphorylation". Albert-Ludwigs-Universität Freiburg im Breisgau, Ph.D. dissertation, p. 72 (March 2017).
37. Putta, Priya, "The Tale/Head of Two Membrane Lipids Through Protein Interactions". Kent State University, Ph.D. dissertation, p. 68 (May 2018).
38. Ansor, Nurhuda Mohamad, "Analysis of the cellular basis for developmental disorders arising from mutations in RAC1". University of Manchester, Ph.D. dissertation, p. 48 (2019).
39. Caygill, Claire Heather, "Investigation of glutamate receptors in human pancreatic  $\beta$ -cells and their role in insulin secretion". University of the West of England, Bristol, Ph.D. dissertation, p. 51 (May 2020).
40. Harlow, E. and Lane, D. *Antibodies: A Laboratory Manual*. Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY), p. 549 (1988).
41. Green, N.M., and Toms, E.J., *Biochem. J.*, **133(4)**, 687-700 (1973).
42. Heney, G., and Orr, G.A., *Anal. Biochem.*, **114(1)**, 92-96 (1981).
43. Orr, G.A., *J. Biol. Chem.*, **256(2)**, 761-766 (1981).

44. Hoffman, W.L., and Jump. A.A., *Anal. Biochem.*, **181(2)**, 318-320 (1989).
45. Houen, G., and Hansen, K., *J. Immunol. Methods*, **210(2)**, 115-123 (1997).
46. Duhamel, R.C., and Whitehead, J.S., *Methods Enzymol.*, **184**, 201-207 (1990).

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

## Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, Millipore, MILLIPLEX and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

S1638pis Rev 01/22 KAT,GCY

