

## Product Information

## Lyticase from *Arthrobacter luteus*


Lyophilized powder,  $\geq 2,000$  units/mg protein, Protein  $\geq 20\%$  by Biuret

**L2524**

### Product Description

CAS Registry Number: 37340-57-1

Yeast cells are difficult to disrupt because the cell walls may form capsules or resistant spores. DNA can be extracted from yeast by using lysing enzymes such as lyticase, chitinase, zymolase, and gluculase to induce partial spheroplast formation. Spheroplasts are subsequently lysed to release DNA. Lyticase is preferred to digest cell walls of yeast and generate spheroplasts from fungi for transformation.<sup>1</sup>



Lyticase contains  $\beta$ -(1 $\rightarrow$ 3)-glucan laminaripentaohydrolase with additional  $\beta$ -(1 $\rightarrow$ 3)-glucanase, protease, and mannanase activities.<sup>1,2</sup> For isolation of nucleic acids, lyticase has been used in the lysis of yeast cell walls (such as from *Candida*, *Debaryomyces*, *Saccharomyces*, *Saccharomycopsis*, *Saccharomycodes*, *Eremothecium*, and *Schwanniomyces* species).<sup>3,4</sup> Several publications cite use of this product in their protocols.<sup>5,6</sup> Several theses<sup>7-11</sup> and dissertations<sup>12-22</sup> have cited use of product L2524 in their protocols.

### Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

### Storage/Stability

Several publications have reported preparation of stock solutions of this product at various concentrations and in different systems, with storage of the stock solutions at  $-20\text{ }^{\circ}\text{C}$ , as follows:

- 20,000 U/mL in ddH<sub>2</sub>O<sup>23</sup>
- 25,000 U/mL in 1 $\times$  PBS<sup>5</sup>

A separate publication reports preparation of stock solutions of lyticase, though not this specific product, at 5 mg/mL in 1 M sorbitol with 0.1 M EDTA, pH 8.0, with storage of this stock solution at  $-20\text{ }^{\circ}\text{C}$  in frozen aliquots.<sup>24</sup> However, we have not tested any of these methods ourselves.

### References

1. Scott, D.H., and Schekman, R., *J. Bacteriol.*, **42(2)**, 414-423 (1980).
2. Machado, F.R.S., Jr. *et al.*, *J. Biotechnol.*, **218**, 108-114 (2016).
3. Van Burik, J.-A.H. *et al.*, *Med. Mycol.*, **36(5)**, 299-303 (1998).
4. Goldschmidt, P. *et al.*, *PLoS One*, **9(6)**, e94886 (2014).
5. Murray, S., and Kiseleva, E., *Methods Cell Biol.*, **88**, 367-387 (2008).
6. Zenklusen, D., and Singer, R.H., *Methods Enzymol.*, **470**, 641-659 (2010).
7. Nelson, Jasen, "Screening of endocrine disrupting compounds using estrogen receptor transcriptional activation *in vitro* bioassays". Simon Fraser University, M.Sc. thesis, p. 84 (2007).
8. Scott, Adam Matthew, "Modified Yeast Two-Hybrid Screening Identifies SKAP-HOM as a Novel Substrate of PTP-PEST". McGill University, M.Sc. thesis, p. 13 (2008).

9. Helminski, Peter Simon, "Understanding the functional role of the *Xenopus laevis* cortical granule lectin oligosaccharides in the block to polyspermic fertilization". California State University, Sacramento, M.S. thesis, p. 23 (2012).
10. Au Yeung, Wan Kin, "Study of the Possible Roles of *OsFKBP12* in Plant Defense System". The Chinese University of Hong Kong, M.Phil. thesis, p. 106 (2013).
11. Chong, Lisa Pui-Ling, "The role of ER-localized HSP90 in modulating the WUSCHEL-CLAVATA pathway and ER stress resistance in *Arabidopsis thaliana*". University of Toronto, M.Sc. thesis, p. 29 (2015).
12. Coque, Laurent, "Function of the *ENOD8* gene in nodules of *Medicago truncatula*". University of North Texas, Ph.D. dissertation, p. 134 (2006).
13. BeGora, Michael D., "Cloning and biochemical characterization of a phosphobase *N*-methyltransferase from *Arabidopsis* involved in choline and phospholipid metabolism". McMaster University, Ph.D. dissertation, p. 58 (2010).
14. Ramsden, Sarah F.L., "The role of gephyrin, collybistin and novel GEFs in the synaptic clustering of inhibitory receptors". University College London, Ph.D. dissertation, p. 54 (2012).
15. Weber, Lisa, "Characterization of *Schizosaccharomyces pombe* Sup11p, a protein involved in  $\beta$ -1,6-glucan biosynthesis". Ruperto-Carola University Heidelberg, Dr. rer. nat. dissertation, p. 24 (2013).
16. Gál, Csenge, "The Impact of the Chromatin Regulators, Abo1 and HIRA, on Global Nucleosome Architecture". Newcastle University, Ph.D. dissertation, p. 70 (2014).
17. Moore, Michael John, "An integrative biophysical and bioanalytical approach for investigating the mitochondrial labile iron pool". Texas A&M University, Ph.D. dissertation, p. 46 (2017).
18. Hampton-Smith, Rachel Jane, "Identification and Characterisation of Novel Substrates and Binding Partners of the Asparaginyl Hydroxylase, FIH". University of Adelaide, Ph.D. dissertation, p. 46 (2015).
19. Sukhareuski, Andrei, "How telomerase and Dna2 govern the fate of chromosome ends in *Saccharomyces cerevisiae*". University of Edinburgh, Ph.D. dissertation, p. 105 (2019).
20. Muraglia, Katrina Anne, "Amphotericin B ion channels restore cystic fibrosis airway surface physiology". University of Illinois at Urbana-Champaign, Ph.D. dissertation, p. 68 (2018).
21. Marr, Luke Timothy, "Probing Genome-Wide Native Chromatin Structure via Histone Surface Accessibilities". University of Rochester, Ph.D. dissertation, pp. 37, 69, 106 (2020).
22. O'Kane, Callum Jude, "Probing the role of histone acetylation in the evolution of pathogenicity in *Candida glabrata*". Queen's University Belfast, Ph.D. dissertation, p. 81 (2020).
23. Lam, F.H. *et al.*, *Protocol Exchange*, doi:10.1038/nprot.2008.121 (2008).
24. Hong, S.-B., Rashid, M.B., and Santiago-Vázquez, L.Z., *Methods in Biotechnology*, "AMB 1 experiment 31: Plasmid DNA isolation from yeast". Wiley-Blackwell (Hoboken, NJ), pp. 143-144 (2017).

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.  
L2524pis Rev 06/22 GCY,DT,LB,MAM

**MERCK**