

## Product Information

## Aminopeptidase from *Aeromonas proteolytica*

Lyophilized powder, 50-150 units/mg protein

**A8200**

### Product Description

CAS Registry Number: 37288-67-8

Enzyme Commission (EC) Number: 3.4.11.10

pI:<sup>1</sup> 3.0-3.5

$\lambda_{\text{max}}$ : 278 nm<sup>1</sup>

Extinction coefficient:  $E^{1\%} = 14.4$  (278.5 nm)<sup>1</sup>

Synonym: AAP

Aminopeptidases are a family of widely distributed proteases which participate in many significant biological processes, such as protein maturation, hormone production, and peptide digestion.<sup>2-4</sup> While several  $\text{Zn}^{2+}$  peptidases are known to contain a single  $\text{Zn}^{2+}$  ion in their active site,<sup>5-7</sup> a few metalloaminopeptidases, including those from bovine lens,<sup>8-11</sup> *Escherichia coli*,<sup>12</sup> *Aeromonas proteolytica*,<sup>13</sup> and *Streptomyces griseus*<sup>14</sup> have been proven by means of X-ray crystallography to contain a dinuclear metal active site.

This AAP product is a metalloenzyme, which contains 2 atoms of  $\text{Zn}^{2+}$  in a single polypeptide with an approximate molecular weight of 29.5 kDa, as determined by sedimentation. Several crystallographic studies of aminopeptidase from *Aeromonas proteolytica* (AAP) have been published.<sup>15-28</sup> Other publications have reported mechanistic studies on AAP.<sup>19-29</sup>

Several dissertations<sup>30-32</sup> have cited use of product A8200 in their research protocols.

### Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

### Preparation Instructions

This product is soluble in water (1 mg/mL) and aqueous buffers.

### Storage/Stability

This enzyme has a high degree of stability, being stable even to temperatures of 70 °C for several hours.<sup>1</sup> Partial inactivation occurs in 8 M urea. Maximum stability and activity are at pH 8.0-8.5.

The enzyme is stable for several years at -20 °C. It may be lyophilized and reconstituted with little loss in activity.

### References

1. Prescott, J.M., and Wilkes, S.H., *Methods Enzymol.*, **45**, 530-543 (1976).
2. Taylor, A., *FASEB J.*, **7(2)**, 290-298 (1993).
3. Taylor, A., *Trends Biochem. Sci.*, **18(5)**, 167-172 (1993).
4. Gonzales, T., and Robert-Baudouy, J., *FEMS Microbiol. Rev.*, **18(4)**, 319-344 (1996).
5. Christianson, D.W., and Lipscomb, W.N., *Acc. Chem. Res.*, **22(2)**, 62-69 (1989).
6. Mangani, S. *et al.*, *Coord. Chem. Rev.*, **120**, 309-324 (1992).
7. Matthews, B.W., *Acc. Chem. Res.*, **21(9)**, 333-340 (1988).
8. Burley, S.K. *et al.*, *Proc. Natl. Acad. Sci. USA*, **87(17)**, 6878 -6882 (1990).
9. Burley, S.K. *et al.*, *J. Mol. Biol.*, **224(1)**, 113-140 (1992).
10. Sträter, N., and Lipscomb, W.N., *Biochemistry*, **34(28)**, 9200 -9210 (1995).

**MERCK**

11. Kim, H., and Lipscomb, W.N., *Adv. Enzymol. Relat. Areas Mol. Biol.*, **68**, 153-213 (1994).
12. Roderick, S.L., and Matthews, B.W., *Biochemistry*, **32(15)**, 3907-3912 (1993).
13. Chevrier, B. *et al.*, *Structure*, **2(4)**, 283-291 (1994).
14. Greenblatt, H.M. *et al.*, *J. Mol. Biol.*, **265(5)**, 620-636 (1997).
15. Chevrier, B. *et al.*, *Structure*, **237(2)**, 393-398 (1996).
16. Stamper, C. *et al.*, *Biochemistry*, **40(24)**, 7035-7045 (2001).
17. Desmarais, W.T. *et al.*, *Structure*, **10(8)**, 1063-1072 (2002).
18. Stamper, C.C. *et al.*, *Biochemistry*, **43(30)**, 9620-9628 (2004).
19. Chen, G. *et al.*, *Biochemistry*, **36(14)**, 4278-4286 (1997).
20. Bennett, B., and Holz, R.C., *Biochemistry*, **36(32)**, 9837-9846 (1997).
21. De Paola, C.C. *et al.*, *Biochemistry*, **38(28)**, 9048-9053 (1999).
22. Ustynyuk, L. *et al.*, *Biochemistry*, **38(35)**, 11433-11439 (1999).
23. Huntington, K.M. *et al.*, *Biochemistry*, **38(47)**, 15587-15596 (1999).
24. Bienvenue, D.L. *et al.*, *J. Inorg. Biochem.*, **78(1)**, 43-54 (2000).
25. Bienvenue, D.L. *et al.*, *J. Biol. Inorg. Chem.*, **7(1-2)**, 129-135 (2002).
26. Bienvenue, D.L. *et al.*, *Biochemistry*, **41(11)**, 3712-3719 (2002).
27. Bennett, B. *et al.*, *J. Am. Chem. Soc.*, **124(44)**, 13025-13034 (2002).
28. Schürer, G. *et al.*, *Biochemistry*, **43(18)**, 5414-5427 (2004).
29. Bzymek, K.P., and Holz, R.C., *J. Biol. Chem.*, **279(30)**, 31018-31025 (2004).
30. Touati, Jeremy, "Phage Selection of Bicyclic Peptide Ligands and Development of a New Peptide Cyclisation Method". École Polytechnique Fédérale de Lausanne, Ph.D. dissertation, p. 68 (2012).
31. Štefková, Kristýna, "Energy metabolism and enzymatic activity in the *Ips typographus* in relation to diapause". University of South Bohemia (České Budějovice), Ph.D. dissertation, p. 55 (2017).
32. Zhao, Xiling, "Discovery and Structural Characterization of Lantipeptides from Soil and Rumen Bacteria". University of Illinois at Urbana-Champaign, Ph.D. dissertation, p. 66 (2017).

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

Merck and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

A8200pis Rev 04/22 CMH,NSB,GCY

