



SIGMA-ALDRICH

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Tetracaine hydrochloride

Product Number **T 7508**
Store at Room Temperature

Product Description

Molecular formula: $C_{15}H_{24}N_2O_2 \cdot HCl$
Molecular weight: 300.8
CAS number: 136-47-0
Melting point: 147-150 °C¹
pK_a: 8.39¹
Extinction coefficients:¹
 E^{mM} (water) = 14.108 (225 nm), 26.352 (310 nm)
 E^{mM} (methanol) = 7.586 (226 nm), 29.512 (310 nm)

Tetracaine hydrochloride is an anesthetic. A study on its mechanism of action² and a comprehensive description³ of this compound have been published. It blocks voltage-sensitive release of Ca^{2+} from sarcoplasmic reticulum. It is an inhibitor of cyclic nucleotide-gated channels.⁴

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

Tetracaine hydrochloride is soluble in water (50mg/ml), yielding a clear, colorless solution. It is also soluble in 40 parts alcohol and 30 parts chloroform. It is practically insoluble in ether, benzene, or acetone.¹

Storage/Stability

Aqueous solutions of tetracaine hydrochloride are stable and can be sterilized by brief boiling.¹ Solutions are incompatible with alkalis. Aqueous solutions should be discarded if they contain crystals or are discolored or cloudy.⁵

References

1. The Merck Index, 12th ed., Entry# 9330.
2. Leung, Y. W., and Rawal, B. D., Mechanism of action of tetracaine hydrochloride against *Pseudomonas aeruginosa*. J. Infect. Dis., **136**(5), 679-683 (1977).
3. Riaz, M., in Analytical Profiles of Drug Substances, Vol. 18, Florey, K., ed., Academic Press (New York, NY: 1989), p. 379-411.
4. www.muhealth.org/~research/srd2001/2001Hanson.PDF
5. Martindale The Extra Pharmacopoeia, 31st ed., Reynolds, J. E. F., ed., Royal Pharmaceutical Society (London, England: 1996), p. 1323.

HLD/RXR 5/06

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.